	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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S4 CHO medium is a chemical defined basal medium with no animal-derived ingredients, and no protein ingredients, which is suitable for batch culture, fed-batch culture and perfusion culture in the development and production of therapeutic protein products by Chinese hamster ovary (CHO). S4 CHO medium is free of hypoxanthine, thymine, L-glutamine and contains small amounts of yeast and soy-derived hydrolysates. Suitable for the culture of CHO-K1 cell lines using GS and DHFR screening systems.

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S4 CHO medium is specially designed to improve the growth and production performance of Chinese hamster ovary cell.

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When cultivating GS screening system engineering cells, hypoxanthine and thymine were added as needed;

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When cultivating DHFR screening system engineering cells, add 2-8mM L-glutamine as needed.

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If insulin is added, it is recommended to add a concentration of 0.5-10mg/L.

【配制指南 Preparation guide】

适用于粉末包装 (以 20L 为例) Suitable for powder package (20L as an example)

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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This medium is low sugar medium, it is recommended to add 6g/L glucose.

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Prepare 120g of glucose, 20ml of 36.5% concentrated hydrochloric acid, 500mL of 0.4M NaOH and 42g of sodium bicarbonate.

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Dissolve 42g of sodium bicarbonate in 1L of purified water.

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Prepare a cleaned 20L glass beaker or similar container and add 10L of purified water at a temperature of about 70°C.

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Take the preparation of 20L S4 CHO dry powder medium - component A (19.82g/L), add it to the water while stirring for 10~15 minutes. When the A powder is completely dissolved, add 1L of solution containing 42g of sodium bicarbonate (dissolve sodium bicarbonate separately before mixing), then add 120g of glucose (corresponding to 6g/L of glucose, which can be adjusted according to the requirements of use) and continue stirring for 10 minutes.

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Take the preparation of 20L S4 CHO dry powder medium - component B (1.48g/L), add 500mL of 0.4M NaOH, stir for 10~15 minutes, and continue to stir for 10 minutes after the B powder is completely dissolved.

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After the temperature of the solution A dropped to room temperature, add prepared 20L of S4 CHO dry powder medium - component C (25mL / L) and completely dissolved B powder solution, stir for 5~10 minutes; add about 20ml of 36.5 % concentrated hydrochloric acid to adjust the pH to 6.9-7.1, and fix the volume to 20L.

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Use a Millipak 40 or similar filter to remove bacteria and filter into a container for use.

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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Most cell lines use this product without any domestication, and can be directly inoculated into this medium and passed for more than three times. For some cell lines, domestication may be used when using this series of medium, and the specific steps are as follows:

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Most cell lines can be directly domesticated to S4 CHO.

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After at least 2~3 generations, the doubling time is normal and stable, and the cell viability is more than 90%, indicating that the cell strain has been domesticated.

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The cell strain was cultured in the original medium to the middle of exponential growth period, and when the cell viability was more than 90%, it was inoculated into the mixed medium with the volume ratio of 50%: 50% (S4 CHO: original medium), and the inoculation density was $3\sim5\times10^5$ cells/mL, and it was cultured at 36.5°C and 6% CO₂. The cells were cultured for 3~4 days to reach more than 1×10^6 cells/mL, and then subcultured.

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Cells were inoculated into a mixed medium with the volume ratio of 75%: 25% (S4 CHO: original medium), and the inoculation density was $3\sim5\times10^5$ cells/mL, and cultured at 36.5°C and 6% CO₂. The cells were cultured for 3~4 days to reach more than 1×10^6 cells/mL, and then subcultured.

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The cells were inoculated into 100% S4 CHO with the inoculation density of $3\sim5 \times 10^5$ cells/mL, and cultured at 36.5°C and 6%CO₂. The cells were cultured for 3~4 days to reach more than 1×10^6 cells/mL, and then subcultured.

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When adopting this domestication procedure, if the cells still grow slowly or have low activity, you can consider slowly increasing the ratio of S4 CHO to 25: 75, 50: 50, 75: 25, 100: 0 from the mixed medium prepared with a volume ratio of 10: 90 (S4 CHO: original medium). Or the cells are collected by centrifugation during the process and subcultured again.

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Prepare frozen solution on the super clean workbench : 90% S4 CHO + 10% dimethyl sulfoxide (DMSO) mixed solution, precooling at 2~8°C (heat will be released when DMSO is diluted);

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Frozen cell fluid: seed cells were in the exponential growth period, the density is greater than 1.5×10^6 cells/mL, and the viability is greater than 95%.

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Cell fluid was centrifuged at 800rpm for 5 min;

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Slowly pour out the supernatant, resuspend the cells with cryopreservation solution, the cryopreservation density is $1.0\sim1.5 \times 10^7$ cells/mL, and transfer the cells to a sterile cryopreservation tube;

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Place the cryopreservation tube in the cryopreservation box containing isopropyl alcohol, freeze it at - 80 °C overnight, and then transfer it to the liquid nitrogen tank for long-term storage. If there is no freezing box, the temperature can be reduced manually by gradient as follows:

- 4°C冻存 30min;
- freeze at 4°C for 30min;
- -20°C冻存 2~4 小时;
- freeze at -20°C for 2~4h;
- -80°C冻存过夜;
- freeze at - 80°C overnight;
- 转移至液氮罐中长期贮存。
- transfer frozen cells to liquid nitrogen tank for long-term storage.

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Prepare a 36.5°C warm water to thaw frozen cells;

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Take out the frozen tube from the liquid nitrogen tank and quickly thaw frozen cells in 36.5°C warm water;

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After wiping the cryopreservation tube with 75% ethanol, open the cryopreservation tube in the sterile operation table, transfer the cell fluid to a 15 ml centrifuge tube containing 2-5 mL of S4 CHO, blow and mix well, centrifuge at 800 rpm for 5 minutes;

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Slowly pour out the supernatant, resuspend with 20~30 ml preheated S4 CHO, and transfer to a 125 ml shake flask;

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Place it in a shaking incubator with 8% CO₂, 80% humidity, 110 ~ 130rpm, at 36.5°C for culture;

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After 2~3 days of culture, the cells were counted and subcultured.

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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The cells are seeded at 5E5 ~6E5, count and subculture every 2 ~ 3 days. In the first three passages, the volume remained unchanged to restore cell viability. After the cell viability recovers to normal and reaches more than 90%. The seed cells were expanded at the density of 5E5 ~6E5 until reaching the required volume. The criteria for normal seed state: the viability was greater than 95%, the cell morphology was regular and round, and the growth doubling time was normal.

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Shanghai production base, S4 CHO dry powder medium - component A, S4 CHO dry powder medium - component B, and S4 CHO dry powder medium - component C: 2°C to 8°C, protect from light; validity period: 12 months.

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Wuxi production base, S4 CHO dry powder medium - component A, and S4 CHO dry powder medium - component B: 2°C to 8°C, protect from light; retest date: 12 months. S4 CHO dry powder medium - component C: 2°C to 8°C, protect from light; validity period: 12 months.

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地址 Address: 上海市松江区民强路 1525 弄 30 幢 Building 30, Lane 1525, Minqiang Road, Songjiang District, Shanghai

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