	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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S3 CHO medium is a chemical defined basal culture medium with no animal-derived ingredients, and no protein ingredients, which is suitable for batch culture, fed-batch culture and perfusion culture in the development and production of therapeutic protein products by Chinese hamster ovary (CHO). S3 CHO medium does not contain phenol red, hypoxanthine, thymine, L-glutamine and contains a small amount of yeast-derived hydrolysate. It is suitable for the culture of different cell lines such as CHO-S and DG44 using GS and DHFR screening systems.

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S3 CHO medium is specially designed to improve the growth and production performance of Chinese hamster ovary cell.

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This product does not contain phenol red, antibiotics, L-glutamine, or insulin;

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When cultivating GS screening system engineering cells, hypoxanthine and thymine were added as needed;

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When cultivating DHFR screening system engineering cells, add 2-8mM L-glutamine as needed.

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If insulin is added, it is recommended to add a concentration of 0.5-10mg/L.

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This medium is low sugar medium, it is recommended to add 6g/L glucose.

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Dissolve 42 g of sodium bicarbonate in 1 L of purified water.

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Prepare a cleaned 20L glass beaker or similar container and add 10L of purified water at a temperature of about 70 °C;

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Take the preparation of 20L of S3 CHO dry powder medium - component A (18.63g/L), add it to the water while stirring for 10~15 minutes. After the A powder is completely dissolved, add 1L containing 42g of sodium bicarbonate solution (dissolve sodium bicarbonate separately before mixing), then add 120g of glucose (corresponding to 6g/L of glucose, which can be adjusted according to the requirements of use), and continue to stir for 10 minutes;

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Take the preparation of 20L of S3 CHO dry powder medium - component B (2.96g/L), add 1L of 0.4M NaOH, stir for 10-15 minutes, and after the B powder is completely dissolved, continue to stir for 10 minutes;

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Take the preparation of 20L of S3 CHO dry powder medium - component C (1.068g/L), add 2L of purified water, stir for 10 minutes, add 9.6mL of 0.4M NaOH, stir for 15 minutes, and obtain the C powder solution;

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After the temperature of the solution A dropped to room temperature, add completely dissolved B and C powder

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

solution, stir for 5~10 minutes; adjust the pH to 6.9-7.1 and fix the volume to 20L.

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Sterilize and filter with a 0.22μm filter.

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Most cell lines use this product without any domestication, and can be directly inoculated into this medium and passed for more than three times. For some cell lines, domestication may be used when using this series of medium, and the specific steps are as follows:

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Most cell lines can be directly domesticated to S3 CHO.

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After at least 2~3 generations, the doubling time is normal and stable, and the cell viability is more than 90%, indicating that the cell strain has been domesticated.

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The cell strain was cultured in the original medium to the middle of exponential growth period, and when the cell viability was more than 90%, it was inoculated into the mixed medium with the volume ratio of 50%: 50% (S3 CHO: original medium), and the inoculation density was  $3\sim5\times10^5$  cells/mL, and it was cultured at 36.5°C and 6% CO<sub>2</sub>. The cells were cultured for 3~4 days to reach more than  $1\times10^6$  cells/mL, and then subcultured.

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Cells were inoculated into a mixed medium with the volume ratio of 75%: 25% (S3 CHO: original medium), and the inoculation density was  $3\sim5\times10^5$  cells/mL, and cultured at 36.5°C and 6% CO<sub>2</sub>. The cells were cultured for 3~4 days to reach more than  $1\times10^6$  cells/mL, and then subcultured.

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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The cells were inoculated into 100% S3 CHO with the inoculation density of  $3\sim5\times10^5$  cells/mL, and cultured at 36.5°C and 6%CO<sub>2</sub>. The cells were cultured for 3~4 days to reach more than  $1\times10^6$  cells/mL, and then subcultured.

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In 100% S3 CHO, at least 2~3 generations, the doubling time is normal and stable, and the cell viability is more than 90%, indicating that the cell strain has been domesticated.

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When adopting this domestication procedure, if the cells still grow slowly or have low activity, you can consider slowly increasing the ratio of S3 CHO to 25: 75, 50: 50, 75: 25, 100: 0 from the mixed medium prepared with a volume ratio of 10: 90 (S3 CHO: original medium). Or the cells are collected by centrifugation during the process and subcultured again.

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Prepare frozen solution on the super clean workbench : 90% S3 CHO + 10% dimethyl sulfoxide (DMSO) mixed solution, precooling at 2~8°C (heat will be released when DMSO is diluted);

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Frozen cell fluid: seed cells were in the exponential growth period, the density is greater than  $1.5\times10^6$  cells/mL, and the viability is greater than 95%.

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Cell fluid was centrifuged at 800rpm for 5 min;

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Slowly pour out the supernatant, resuspend the cells with cryopreservation solution, the cryopreservation density is  $1.0\sim1.5\times10^7$  cells/mL, and transfer the cells to a sterile cryopreservation tube;

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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Place the cryopreservation tube in the cryopreservation box containing isopropyl alcohol, freeze it at - 80 °C overnight, and then transfer it to the liquid nitrogen tank for long-term storage. If there is no freezing box, the temperature can be reduced manually by gradient as follows:

- 4℃冻存 30min;
- freeze at 4℃ for 30min;
- -20℃冻存 2~4 小时;
- freeze at -20℃ for 2~4h;
- -80℃冻存过夜;
- freeze at - 80℃ overnight;
- 转移至液氮罐中长期贮存。
- transfer frozen cells to liquid nitrogen tank for long-term storage.

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Prepare a 37.0℃ warm water to thaw frozen cells;

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Take out the frozen tube from the liquid nitrogen tank and quickly thaw frozen cells in 36.5℃ warm water;

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After wiping the cryopreservation tube with 75% ethanol, open the cryopreservation tube in the sterile operation table, transfer the cell fluid to a 15 ml centrifuge tube containing 2-5 mL of S3 CHO, blow and mix well, centrifuge at 800 rpm for 5 minutes;

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Slowly pour out the supernatant, resuspend with 20~30 ml preheated S3 CHO, and transfer to a 125 ml shake flask;

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Place it in a shaking incubator with 8% CO<sub>2</sub>, 80% humidity, 110 ~ 130rpm, at 36.5℃ for culture;

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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After 2~3 days of culture, the cells were counted and subcultured.

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The cells are seeded at 5E5 ~6E5, count and subculture every 2 ~ 3 days. In the first three passages, the volume remained unchanged to restore cell viability. After the cell viability recovers to normal and reaches more than 90%. The seed cells were expanded at the density of 5E5 ~6E5 until reaching the required volume. The criteria for normal seed state: the viability was greater than 95%, the cell morphology was regular and round, and the growth doubling time was normal.

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Wuxi production base, S3 CHO dry powder medium - component A, S3 CHO dry powder medium - component B, and S3 CHO dry powder medium - component C: 2°C to 8°C, protect from light; retest date: 12 months.

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