	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

## 产品使用说明书 Product Use Instructions

多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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2×10<sup>7</sup> is a kind of high-efficiency fed-batch culture medium without animal-derived components and protein, which is suitable for the development and production of therapeutic protein products by Chinese hamster ovary (CHO). This product contains no hypoxanthine, thymine, L- glutamine, P188 and sodium bicarbonate. It is suitable for the culture of different cell lines, such as CHO-K1, CHO-S, DG44 and CHOZN, using GS system and dhfr screening system. 2×10<sup>7</sup> is divided into 2×10<sup>7</sup> high-efficiency feed medium A and high-efficiency feed medium B.

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2×10<sup>7</sup> is an efficient feeding formula specially designed to improve the growth and production performance of CHO cells. It is rich in nutrients such as amino acids and vitamins, which can meet the nutritional needs of high-density culture of cells and high expression of products. Combined with the characteristics of cell lines, reasonable feeding strategy can obviously improve the growth of cells and antibody expression.

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It is suggested that the first feeding should be started in the middle and late period of cell exponential growth, when

the cell density is above  $2.0 \times 10^6$  cells/mL. The recommended feeding amount is: 5% on the 3rd day, 5% on the 5th day, 7.5% on the 7th day, 7.5% on the 9th day, 5% on the 11th day and 5% on the 13th day.

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It is suggested that the total feeding amount should not exceed 40%.

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When cultivating GS system cell lines, hypoxanthine and thymine were added as needed.

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When the cell lines of dhfr system are cultured, 3~8mM L- glutamine is added as needed.

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The liquid packaging of this product will produce some precipitation when it is placed at 2~8 °C , so it is recommended to use it within two months.

### 【配制指南 Preparation guide】

适用于粉末包装（以 1L 为例）Suitable for powder packaging (taking 1L as an example)

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Add 107.09g of  $2 \times 10^7$  high-efficiency feed medium A powder, stir for 20 min, and dissolve completely. At this time, the pH is around 4.0;

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Remarks: If the powder of  $2 \times 10^7$  high-efficiency feeding medium A is not clear after mixing, the temperature can be controlled in the process, or the temperature can be appropriately increased (the maximum is not more than 40°C) until the solution is completely clear.

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In addition, prepare 50 ml of 0.4 M NaOH solution, add 2.96g of high-efficiency feed medium B powder, stir for 5 min, and fully dissolve it, and record it as solution B;

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Add solution B and 1g P188 in turn, and continue stirring for 10min;

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The normal osmotic pressure of 2 $\times$ 10<sup>7</sup> medium is 800~900 mOsm/kg.

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Shanghai production base, powder packaging: 2 $^{\circ}$ C to 8 $^{\circ}$ C, protect from light; validity period: 24 months.

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Wuxi production base, powder packaging: 2 $^{\circ}$ C to 8 $^{\circ}$ C, protect from light; retest date: 24 months.

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地址 Address: 上海市松江区民强路 1525 弄 30 幢 Building 30, Lane 1525, Minqiang Road, Songjiang District, Shanghai

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