	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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2×103A high-efficiency feed medium is a kind of high-efficiency fed-batch culture medium without animal-derived components and protein, which is suitable for the development and production of therapeutic protein products by Chinese hamster ovary (CHO). This product contains no hypoxanthine, thymine, L- glutamine, P188 and sodium bicarbonate. It is suitable for the culture of different cell lines, such as CHO-K1, CHO-S, DG44 and CHOZN, using GS system and dhfr screening system. 2×103A high-efficiency feed medium is divided into 2×103A high-efficiency feed medium A and high-efficiency feeding medium B.

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2×103A high-efficiency feed medium is an efficient feeding formula specially designed to improve the growth and production performance of CHO cells. It is rich in nutrients such as amino acids and vitamins, which can meet the nutritional needs of high-density culture of cells and high expression of products. Combined with the characteristics of cell lines, reasonable feeding strategy can obviously improve the growth of cells and antibody expression.

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It is suggested that the first feeding should be started in the middle and late period of cell exponential growth, when

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

the cell density is above  $2.0 \times 10^6$  cells/mL. The recommended feeding amount is: 5% on the 3rd day, 5% on the 5th day, 7.5% on the 7th day, 7.5% on the 9th day, 5% on the 11th day and 5% on the 13th day.

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It is suggested that the total feeding amount should not exceed 40%.

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2×103A high-efficiency feed medium can be used in combination with BR7, and the culture efficiency can be further improved. The recommended concentration of BR7 for each addition is 1%.

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Hypoxanthine and thymine were added as needed when cultivating GS screening system engineering cell lines.

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3~8mM L- glutamine was added as needed when cultivating dhfr screening system engineering cell lines.

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The liquid packaging of this product will produce some precipitation when it is placed at 2~8 °C , so it is recommended to use it within two months.

#### 【配制指南 Preparation guide】

适用于粉末包装（以 1L 为例） Suitable for powder packaging (taking 1L as an example)

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Add 99.03g of 2×103A high-efficiency feed medium A powder, stir for 20 min, and dissolve completely. At this time, the solution is slightly turbid, and the pH is about 4.0;

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In addition, prepare 50 ml of 0.4M NaOH solution, add 2.96g of high-efficiency feeding medium B powder, stir for 10 min, and fully dissolve it, and record it as solution B;

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Add solution B and continue stirring for 20 min;

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After adjusting pH to about 6.70, the solution was clarified.

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The normal osmotic pressure of 2 $\times$ 103A high-efficiency feed culture medium is 750~850 mOsm/kg.

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Shanghai production base, powder packaging: 2 $^{\circ}$ C to 8 $^{\circ}$ C, protect from light; validity period: 24 months.

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Wuxi production base, powder packaging: 2 $^{\circ}$ C to 8 $^{\circ}$ C, protect from light; retest date: 24 months.

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地址 Address: 上海市松江区民强路 1525 弄 30 幢 Building 30, Lane 1525, Minqiang Road, Songjiang District, Shanghai

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