	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNING

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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BR7 CHO super concentrated high-efficiency feed medium is a chemical defined high-efficiency feeding medium without animal-derived ingredients or protein, which is suitable for fed-batch culture in the research and development, and production of therapeutic protein products using Chinese hamster ovary (CHO).

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BR7 CHO needs to be used in combination with AR7 CHO super concentrated high-efficiency feed medium, and can also be used in combination with DuoNING feed medium (such as 2×103A high-efficiency feed medium) and other commercial feed medium, which can further improve the cell growth and expression.

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BR7 CHO super concentrated high-efficiency feed medium is an efficient formula specially designed to further improve the growth and production performance of CHO cells. Its pH is alkaline and does not contain glucose. It should be used in combination with AR7 CHO super concentrated high-efficiency feed medium or other feed medium (such as 2×103A high-efficiency feed medium), and the recommended addition ratio is 1%.

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Suitable for BR7 CHO super concentrated high-efficiency feed medium (powder packaging: taking 1L as an example)

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Add 19.85g of BR7 CHO super concentrated high-efficiency feed medium powder, stir for 30 min, and dissolve completely;

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Sterilize and filter with a 0.22μm filter;

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The osmotic pressure range of BR7 CHO super-concentrated high-efficiency feed liquid is 250~350 mOsm/kg.

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Shanghai production base, powder packaging: 2°C to 8°C, protect from light; validity period: 24 months.

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Wuxi production base, powder packaging: 2°C to 8°C, protect from light; retest date: 24 months.

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地址 Address: 上海市松江区民强路 1525 弄 30 幢 Building 30, Lane 1525, Minqiang Road, Songjiang District, Shanghai

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