记录编号 FN: SMP-RD-A0002-R06

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归档 Filing: filed immediately, permanently kept

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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Media C-03 is a chemical defined basal medium with no animal origin components, no protein, which is suitable for batch culture, fed-batch culture and perfusion culture in the development and production of therapeutic protein products by Chinese hamster ovary (CHO). Media C-03 does not contain L- glutamine. It is suitable for the culture of different cell lines, such as CHOK1, CHOK1SV, CHOS, DG44.

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Most cell lines use this product without any domestication, and can be directly inoculated into this medium and passed for more than three times. For some cell lines, domestication may be used when using this series of medium, and the specific steps are as follows:

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Most cell lines can be directly domesticated to Media C-03.

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After at least 2~3 generations, the doubling time is normal and stable, and the cell viability is more than 90%, indicating that the cell strain has been domesticated.

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The cell strain was cultured in the original medium to the middle of exponential growth period, and when the cell viability was more than 90%, it was inoculated into the mixed medium with the volume ratio of 50%: 50% (Media C-03: original medium), and the inoculation density was $3\sim 5\times 10^5$ cells/mL, and it was cultured at 36.5°C and 6% CO₂. The cells were cultured for 3~4 days to reach more than 1×10^6 cells/mL, and then subcultured.

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Cells were inoculated into a mixed medium with the volume ratio of 75%: 25% (Media C-03: original medium), and the inoculation density was $3\sim 5\times 10^5$ cells/mL, and cultured at 36.5°C and 6% CO₂. The cells were cultured for 3~4 days to reach more than 1×10^6 cells/mL, and then subcultured.

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The cells were inoculated into 100% Media C-03 with the inoculation density of $3\sim 5\times 10^5$ cells/mL, and cultured at 36.5°C and 6%CO₂. The cells were cultured for 3~4 days to reach more than 1×10^6 cells/mL, and then subcultured.

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In 100% Media C-03, at least 2~3 generations, the doubling time is normal and stable, and the cell viability is more than 90%, indicating that the cell strain has been domesticated.

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When adopting this domestication procedure, if the cells still grow slowly or have low activity, we can consider slowly increasing the ratio of Media C-03 to 25: 75, 50: 50, 75: 25, 100: 0 from the mixed medium prepared with

a volume ratio of 10: 90 (Media C-03: original medium). Or the cells are collected by centrifugation during the process and subcultured again.

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Prepare frozen solution on the super clean workbench : 90% Media C-03 +10% dimethyl sulfoxide (DMSO) mixed solution, precooling at 2~8°C (heat will be released when DMSO is diluted);

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Frozen cell fluid: in the exponential growth period, the density is greater than 1.5×10^6 cells/mL, and the viability is greater than 95%. Generally, it is recommended that the frozen storage density is $1.0 \sim 1.5 \times 10^7$ cells/ml;

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Cell fluid was centrifuged at 800rpm for 5 min;

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Slowly pour out the supernatant, resuspend the cells with cryopreservation solution, the cryopreservation density is $1.0 \sim 1.5 \times 10^7$ cells/mL, and transfer the cells to a sterile cryopreservation tube;

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Place the cryopreservation tube in the cryopreservation box containing isopropyl alcohol, freeze it at - 80 °C overnight, and then transfer it to the liquid nitrogen tank for long-term storage. If there is no freezing box, the temperature can be reduced manually by gradient as follows:

- 4°C 冻存 30min;
- freeze at 4°C for 30min;
- -20°C 冻存 2~4 小时;
- freeze at -20°C for 2~4h;
- -80°C 冻存过夜;
- freeze at - 80°C overnight;

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

- 转移至液氮罐中长期贮存。
- transfer frozen cells to liquid nitrogen tank for long-term storage.

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Prepare a 36.5 °C warm water to thaw frozen cells;

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Take out the frozen tube from the liquid nitrogen tank and quickly thaw the cells in 36.5°C warm water;

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After wiping the cryopreservation tube with 75% ethanol, open the cryopreservation tube in the sterile operation table, transfer the cell fluid to a 15 ml centrifuge tube containing 2-5 mL of Media C-03, blow and mix well, centrifuge at 800 rpm for 5 minutes;

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Slowly pour out the supernatant, resuspend with 20~30 ml preheated Media C-03, and transfer to a 125 ml shake flask;

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Place it in a shaking incubator with 6~8% CO₂, 110 ~ 130rpm, at 36.5°C for culture;

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After 2~3 days of culture, the cells were counted and subcultured.

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The cells are seeded at 5E5 ~6E5, count and subculture every 2 ~ 3 days. In the first three passages, the volume remained unchanged to restore cell viability. After the cell viability recovers to normal and reaches more than 90%. The seed cells were expanded at the density of 5E5 ~6E5 until reaching the required volume. The criteria for normal seed state: the viability was greater than 95%, the cell morphology was regular and round, and the growth doubling time was normal.

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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Media C-03 liquid packaging: 2°C to 8°C, protect from light; Shelf life: 12 months.

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